


## Flagship™


### Frontline Coronary Guidewire

The Flagship interventional coronary guidewire is designed for enhanced device delivery during Percutaneous Transluminal Coronary Angioplasty (PTCA).



- **Proprietary Twister® core technology** to provide 1:1 torque control for optimal steerability
- **Hydrophilic-coated distal tip** to facilitate superior wire trackability and lesion crossing
- **3 cm tip radiopacity** to enable fluoroscopic visualization of the distal tip for precise guidewire positioning

## Product Details



## Ordering Information

Flagship™ Frontline Coronary Guidewire				
Part Number	Diameter (in)	Length (cm)	Support	Tip Shape
1603-000-010	0.014	195*	Normal	Straight (Shapeable)
1603-010-010	0.014	300	Normal	Straight (Shapeable)
1603-001-010	0.014	195*	Normal	J-Tip
1603-011-010	0.014	300	Normal	J-Tip
1603-100-010	0.014	195*	Extra Support	Straight (Shapeable)
1603-110-010	0.014	300	Extra Support	Straight (Shapeable)
1603-101-010	0.014	195*	Extra Support	J-Tip
1603-111-010	0.014	300	Extra Support	J-Tip

\*195 cm length guidewires are extendible.

## Regulatory Support

The Flagship Guidewire line is 510(k) and CE mark registered. All facilities producing the Flagship Guidewire line are ISO 13485 compliant. Integer can provide technical dossier assistance to support customer regulatory submissions.

## Indications for Use

The Flagship guidewire is intended for use in the coronary and peripheral vasculature.



### Partner with the Industry Guidewire Expert

As an established leader, we have over 40 years of guidewire manufacturing experience and produce over 30 million guidewires annually. We are here to help you reduce development costs and expedite time to market through our deep portfolio of market-ready guidewires and custom development solutions. E-mail us at [info@synlas.de](mailto:info@synlas.de) to connect with one of our representatives.