



# Introducer Sheaths



Greatbatch Medical™

# RadialSeal™ Introducer Kit\*

for Transradial Access



Guidewire Compatibility (in)	French Size Range (Fr)	Effective Length (cm)	Configurations
.021	4, 5, 6	10, 16, 25	Single-wall puncture, double-wall puncture

**Single-Wall Puncture Kit:**

- Introducer
- .021 in stainless steel wire (45 & 80 cm length)
- 21-gauge needle (4 cm length)

**Double-Wall Puncture Kit:**


- Introducer
- .021 in plastic guidewire (45 & 80 cm length)
- 22-gauge IV catheter (2.5 cm length)

Private label as well as custom and standard kit configurations available upon request.



Greatbatch Medical™

# OptiSeal® Valved PTFE Peelable Introducer




Guidewire Compatibility (in)	French Size Range (Fr)	Effective Length (cm)	Configurations
.035	6, 7, 8, 9, 10, 11, 12	13, 25	Bulk (Non-sterile), Sterile Kits


**Sterile Kits include:**

- Introducer
- 0.035 in guidewire
- 18 gauge needle
- Syringe

Available in sideport and non-sideport options.  
Private label as well as custom and standard kit configurations available upon request.



Greatbatch Medical™  
**ViaPeel®**  
 PTFE Peelable Introducer




Guidewire Compatibility (in)	French Size Range (Fr)	Effective Length (cm)	Configurations
.018	3.5, 4.5, 5, 6	5, 7, 10	Bulk (Non-Sterile)

ViaPeel is a non-valved introducer  
 Private label as well as custom and standard kit configurations available upon request.

Greatbatch Medical™

# FlowGuard™

## Valved PTFE Peelable Introducer




Guidewire Compatibility (in)	French Size Range (Fr)	Effective Length (cm)	Configurations
.038, .035	7	13	Bulk (Non-Sterile)
.038, .035	8, 9, 10, 13	15	Bulk (Non-Sterile)
.038, .035	14, 15, 16, 17	14	Bulk (Non-Sterile)

Private label as well as custom and standard kit configurations available upon request.

Greatbatch Medical™

# PTFE Peelable Introducer




Guidewire Compatibility (in)	French Size Range (Fr)	Effective Length (cm)	Configurations
.038, .035	7, 8, 9, 10, 12, 13, 14	14	Bulk (Non-Sterile)
.038, .035	15, 16	13	Bulk (Non-Sterile)
.018	3.5, 4.5, 5, 6, 7	5, 10	Bulk (Non-Sterile)

**Sterile Kits include:**

- Introducer
- 0.035 in access wire
- 18 gauge needle
- Syringe

Private label as well as custom and standard kit configurations available upon request.



Greatbatch Medical™

# Coaxial Micro Introducer



Guidewire Compatibility (in)	French Size Range (Fr)	Effective Length (cm)	Configurations
.018	4, 5	9	Bulk (Non-Sterile), Sterile Kits, Individual Pouch (Sterile)
.018	4.5	9	Bulk (Non-Sterile)

**Sterile Kits include:**


- Introducer
- 0.018 in access wire
- 21 gauge needle

Private label as well as custom and standard kit configurations available upon request.

Greatbatch Medical™

# CoAxxess™ Plus

## Stiffer Coaxial Micro Introducer



Guidewire Compatibility (in)	French Size Range (Fr)	Effective Length (cm)	Configurations
.018	4, 5	9	Bulk (Non-Sterile)

**Sterile Kits include:**

- Introducer
- 0.018 in access wire
- 21 gauge needle

Private label as well as custom and standard kit configurations available upon request.



