


# Vulcan LT and ST

Provides excellent crossing performance




## Key Product Benefits

- Excellent Tracking & Torque Response for Vulcan LT (12 gm)
- Strong crossing with Vulcan ST (30 gm)
- High Tactile Feel
- Excellent Crossing

# Vulcan ST design specifications


## Vulcan ST 30 gm (0.014" O.D.)

195 / 300 cm



## Vulcan ST 30 gm (0.018" O.D.)

195 / 300 cm



## Vulcan P/N



P/N	Description	Diameter	Length	Tip Load
901023-05	VULCAN LT, 12G 195-014	0.014"	195cm	12g
901023-06	VULCAN LT, 12G 300-014	0.014"	300cm	12g
901023-07	VULCAN LT, 12G 195-018	0.018"	195cm	12g
901023-08	VULCAN LT, 12G 300-018	0.018"	300cm	12g
901023-09	VULCAN ST, 30G 195-014	0.014"	195cm	30g
901023-10	VULCAN ST, 30G 300-014	0.014"	300cm	30g
901023-11	VULCAN ST, 30G 195-018	0.018"	195cm	30g
901023-12	VULCAN ST, 30G 300-018	0.018"	300cm	30g